

ANALYSIS OF FLUORIDE IN WATER

A rapid and simple method for the measurement of fluoride concentration in water in the range of 0.1 - 100 ppm.

Equipment Required

1. EDT Model DR359TX Ion Meter or pH meter with millivoly Scale
2. EDT directION Fluoride Combination ion selective electrode(3221)

Reagents

1. 100 ppm fluoride standard
2. Total ionic strength adjustment buffer (TISAB)
3. 1 mol/l KNO_3 reference electrode filling solution

Sample Preparation

To a volume of sample add an equal volume of TISAB.

Standard Preparation

Make up standards of 100, 10, 1 and 0.1 ppm by serial dilution of the 100 ppm standard with deionised water.

Dilute each standard with an equal volume of TISAB.

Method

Immerse the electrode tip in each of the standards, least concentrated first and most concentrated last, rinsing the electrodes with distilled water between standards.

Plot a graph on lin/log graph paper of mV response against standard concentration.

Immerse the electrode tip in the sample and plot sample concentration from the graph.

This determination may be carried out directly in concentration units by use of the "concentration" modes on EDT directION ION Meters

Calculations

As both standards and samples have been diluted by the same amount the result obtained from the graph is the concentration of the original sample.