



## **Flexible Stand/Electrode Holder - E8060**

---

The E8060 is a strong flexible electrode stand with a sprung arm and a cast base for stability.

The holder suits all standard electrodes up to 12mm diameter bodies with up to 16mm diameter caps and includes a hole for the ATC/Temperature probe as required.

The stand is particularly useful when using Glass pH and reference electrodes as this can prevent the sensitive glass making an impact on the bottom of the beaker. It is also essential when using a magnetic stirrer and flea to stir the standards and samples.

For monitoring the pH of a solution in large containers the flexible arm has an extended reach enabling use in 5 and 10 litre beakers/containers for unattended use.

In some instances pH electrodes will suffer static effects from handheld operation or where electrode screening is insufficient. Using the E8060 will eliminate these static effects.

The E8060 Flexible stand/arm is compatible with EDT directlON standard electrodes and most alternative makes.

## SPECIFICATION TABLE

For more information on this product visit [www.edt.co.uk/E8060](http://www.edt.co.uk/E8060)

## Related Products

---



Basic pH Meter



Basic Conductivity Meter



Glass Combination pH Electrode

## Stay in touch with EDT

---



Telephone +44 (0)1304 829960  
The Old Silo Store, St. Radigund's Abbey CT15 7DL UK  
[www.edt.co.uk](http://www.edt.co.uk)

